

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.12	
Mean	177	Mean	263.48
SD	2.77	SD	2.95
CV	1.56	CV	1.12
Day-Run	Value	Day-Run	Value
1 - 1	180	1 - 1	260
2 - 1	179	2 - 1	267
4 - 1	177	4 - 1	263
6 - 1	180	6 - 1	267
7 - 1	177	7 - 1	264
8 - 1	174	8 - 1	267
11 - 1	173	11 - 1	260
12 - 1	177	12 - 1	265.1
13 - 1	179	13 - 1	263
14 - 1	173	14 - 1	260
18 - 1	178	18 - 1	262
19 - 1	180	19 - 1	268
20 - 1	177	20 - 1	260
21 - 1	180	21 - 1	260
22 - 1	180	22 - 1	265
26 - 1	178	26 - 1	266
27 - 1	176	27 - 1	265
28 - 1	172	28 - 1	259
29 - 1	173	29 - 1	265

Westgard	تفسیر	Cholesterol	آزمایش	اتوآنالایزر BS-800	دستگاه	6	ماه	1402	سال
Normal Control	biorex	Mean	176.5	SD	2.41	Exp. Date	1403/1/17	Lot. No	035

