

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	127	
2 - 1	127	
4 - 1	126	
5 - 1	127	
6 - 1	124	
7 - 1	129	
8 - 1	125	
11 - 1	124	
12 - 1	127	
13 - 1	127	
14 - 1	124	
18 - 1	126	
19 - 1	128	
20 - 1	124	
21 - 1	125	
22 - 1	129	
26 - 1	128	
27 - 1	127	
28 - 1	127.2	
29 - 1	125	
29 - 1	125	

Westgard	تفسیر	Copper (Cu)	آزمایش	اتوآنالایزر BS-800	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	124.88	SD	2.89	Exp. Date	1402/9/20	Lot. No	hsn.35
Next Month	Mean	126.25	SD	1.61	CV	1.28			

