

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	2.57
SD	0.05
CV	2.11
Day-Run	Value
1 - 1	2.48
2 - 1	2.5
4 - 1	2.58
5 - 1	2.62
6 - 1	2.5
7 - 1	2.6
8 - 1	2.62
11 - 1	2.6
12 - 1	2.55
13 - 1	2.5
14 - 1	2.55
18 - 1	2.6
19 - 1	2.57
20 - 1	2.64
21 - 1	2.55
22 - 1	2.7
26 - 1	2.56
27 - 1	2.56
28 - 1	2.61
29 - 1	2.55

High - CV : 1.93	
Mean	3.65
SD	0.07
CV	1.93
Day-Run	Value
1 - 1	3.5
2 - 1	3.61
4 - 1	3.7
5 - 1	3.7
6 - 1	3.6
7 - 1	3.71
8 - 1	3.7
11 - 1	3.65
12 - 1	3.7
13 - 1	3.74
14 - 1	3.5
18 - 1	3.7
19 - 1	3.57
20 - 1	3.62
21 - 1	3.7
22 - 1	3.72
26 - 1	3.61
27 - 1	3.65
28 - 1	3.7
29 - 1	3.6

Westgard	تفسیر	Creatinine	آزمایش	اتوآنالایزر BS-800	دستگاه	6	ماه	1402	سال
Normal Control	biorex	Mean	2.5	SD	0.1	Exp. Date	1403/1/17	Lot. No	035

