

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 3.9	
Mean	0.46	Mean	2.18
SD	0.02	SD	0.08
CV	5.4	CV	3.9
Day-Run	Value	Day-Run	Value
3 - 1	0.46	3 - 1	2.15
4 - 1	0.43	4 - 1	2.3
5 - 1	0.47	5 - 1	2.1
6 - 1	0.42	6 - 1	2.2
8 - 1	0.46	8 - 1	2.23
9 - 1	0.47	9 - 1	2.1
10 - 1	0.46	10 - 1	2.1
12 - 1	0.48	12 - 1	2.18
13 - 1	0.51	13 - 1	2.31
15 - 1	0.46	15 - 1	2.17
16 - 1	0.45	16 - 1	2.18
17 - 1	0.46	17 - 1	2.1
20 - 1	0.42	20 - 1	2.3
22 - 1	0.42	22 - 1	2.16
23 - 1	0.46	23 - 1	2.1
24 - 1	0.41	24 - 1	2.32
25 - 1	0.49	25 - 1	2.05
26 - 1	0.45	26 - 1	2.21
27 - 1	0.46	27 - 1	2.17
29 - 1	0.46	29 - 1	2.23
30 - 1	0.47	30 - 1	2.03

Westgard	تفسیر	D.Bili	آزمایش	اتوآنالایزر BS-800	دستگاه	7	ماه	1402	سال
Normal Control	biorex	Mean	0.46	SD	0.03	Exp. Date	1403/1/17	Lot. No	035

