

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	0.46
SD	0.03
CV	6.61
Day-Run	Value
1 - 1	0.46
2 - 1	0.48
4 - 1	0.48
5 - 1	0.41
6 - 1	0.48
7 - 1	0.5
8 - 1	0.45
8 - 1	0.45
11 - 1	0.53
12 - 1	0.47
13 - 1	0.43
14 - 1	0.47
18 - 1	0.46
19 - 1	0.47
20 - 1	0.42
21 - 1	0.5
22 - 1	0.42
27 - 1	0.48
28 - 1	0.43
29 - 1	0.45

High - CV : 4.02	
Mean	2.16
SD	0.09
CV	4.02
Day-Run	Value
1 - 1	2.2
2 - 1	2.1
4 - 1	2.32
5 - 1	2.11
6 - 1	2.18
7 - 1	2.1
8 - 1	2.28
11 - 1	2.05
12 - 1	2.15
13 - 1	2.2
14 - 1	2
18 - 1	2.15
19 - 1	2.2
20 - 1	2.15
21 - 1	2.12
22 - 1	2.17
26 - 1	2.15
27 - 1	2.28
28 - 1	2
29 - 1	2.25

Westgard	تفسیر	D.Bili	آزمایش	اتوآنالایزر BS-800	دستگاه	6	ماه	1402	سال
Normal Control	biorex	Mean	0.46	SD	0.05	Exp. Date	1403/1/17	Lot. No	035

