

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 1.38	
Mean	73.33	Mean	230.14
SD	1.11	SD	3.17
CV	1.51	CV	1.38
Day-Run	Value	Day-Run	Value
3 - 1	72	3 - 1	227
4 - 1	74	4 - 1	230
5 - 1	74	5 - 1	228
6 - 1	74	6 - 1	227
8 - 1	74	8 - 1	233
9 - 1	75	9 - 1	225
10 - 1	72	10 - 1	230
12 - 1	75	12 - 1	229
13 - 1	73	13 - 1	231
15 - 1	72	15 - 1	232
16 - 1	72	16 - 1	235
17 - 1	73	17 - 1	235
20 - 1	73	20 - 1	231
22 - 1	74	22 - 1	230
23 - 1	75	23 - 1	233
24 - 1	75	24 - 1	229
25 - 1	73	25 - 1	228
26 - 1	73	26 - 1	227
27 - 1	72	27 - 1	232
29 - 1	73	29 - 1	225
30 - 1	72	30 - 1	236

Westgard	تفسیر	FBS	آزمایش	BS-800 اتوآنالایزر	دستگاه	7	ماه	1402	سال
Normal Control	biorex	Mean	73.45	SD	1.1	Exp. Date	1403/1/17	Lot. No	035

