

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.62	
Mean	73.45	Mean	229.8
SD	1.1	SD	3.72
CV	1.5	CV	1.62
Day-Run	Value	Day-Run	Value
1 - 1	74	1 - 1	227
2 - 1	75	2 - 1	233
4 - 1	75	4 - 1	232
5 - 1	74	5 - 1	235
6 - 1	73	6 - 1	228
7 - 1	74	7 - 1	232
8 - 1	75	8 - 1	226
11 - 1	72	11 - 1	232
12 - 1	73	12 - 1	230
13 - 1	72	13 - 1	232
14 - 1	73	14 - 1	225
18 - 1	74	18 - 1	234
19 - 1	74	19 - 1	230
20 - 1	73	20 - 1	234
21 - 1	72	21 - 1	234
22 - 1	74	22 - 1	227
26 - 1	72	26 - 1	231
27 - 1	75	27 - 1	228
28 - 1	73	28 - 1	223
29 - 1	72	29 - 1	223

Westgard	تفسیر	FBS	آزمایش	اتوآنالایزر BS-800	دستگاه	6	ماه	1402	سال
Normal Control	biorex	Mean	73.73	SD	0.83	Exp. Date	1403/1/17	Lot. No	035

