

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	25.7	
3 - 1	24.95	
4 - 1	23.71	
5 - 1	24.15	
6 - 1	25.43	
7 - 1	24.94	
8 - 1	25.12	
9 - 1	26.4	
10 - 1	25.9	
12 - 1	25.1	
13 - 1	24.9	
15 - 1	25.9	
16 - 1	26.8	
17 - 1	26.5	
18 - 1	26.1	
19 - 1	25.4	
24 - 1	23.4	
25 - 1	23.8	
27 - 1	25.6	
29 - 1	21.6	1:3S
30 - 1	23.6	

Westgard	تفسیر	FSH	آزمایش	2000 ایمولایت XPI-	دستگاه	7	ماه	1402	سال
Control	Seronorm	Mean	25.19	SD	0.9	Exp. Date	1403/2/10	Lot. No	2011909
Next Month	Mean	25.17	SD	1.01	CV	4			

