

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	157	
2 - 1	158.1	
4 - 1	159.5	
5 - 1	159.5	
6 - 1	164.2	
7 - 1	167.1	
8 - 1	158.4	
9 - 1	164.8	
11 - 1	167.2	
12 - 1	159.4	
13 - 1	156.2	
14 - 1	149.5	
16 - 1	167.9	
18 - 1	165.5	
19 - 1	147	
20 - 1	149.8	
21 - 1	167	
22 - 1	169	
26 - 1	154	
27 - 1	159	
28 - 1	162	
29 - 1	161	
30 - 1	172	

Westgard	تفسیر	Ferritin	آزمایش	2000 ایمولایت XPI-	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	165.5	SD	12.99	Exp. Date	1403/2/10	Lot. No	365-02-028
Next Month	Mean	160.66	SD	6.59	CV	4.1			

