

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	7.3	
3 - 1	7.52	
4 - 1	7.42	1:2S
5 - 1	7.32	
6 - 1	7.31	Cusum 3S
8 - 1	7.45	
9 - 1	7.49	
10 - 1	7.3	
12 - 1	7.43	
13 - 1	7.51	
15 - 1	7.47	
16 - 1	7.5	
17 - 1	7.35	
18 - 1	7.29	
19 - 1	7.35	
20 - 1	7.41	
22 - 1	7.46	
23 - 1	7.53	
24 - 1	7.48	
25 - 1	7.39	
26 - 1	7.35	
27 - 1	7.28	
29 - 1	7.29	

Westgard	تفسیر	Free T3	آزمایش	SIEMENS Centaur CP	دستگاه	7	ماه	1402	سال
Control	biorex	Mean	7.44	SD	0.1	Exp. Date	1401/12/25	Lot. No	36-02-026
Next Month	Mean	7.4	SD	0.09	CV	1.16			

