

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	7.37	
2 - 1	7.29	
4 - 1	7.37	Cusum 3S
5 - 1	7.49	
7 - 1	7.44	
8 - 1	7.52	
9 - 1	7.38	
11 - 1	7.43	
12 - 1	7.47	
13 - 1	7.53	
14 - 1	7.65	
16 - 1	7.45	
18 - 1	7.28	
19 - 1	7.47	
20 - 1	7.39	
21 - 1	7.65	
22 - 1	7.52	
26 - 1	7.35	
27 - 1	7.56	
28 - 1	7.32	
29 - 1	7.42	
30 - 1	7.38	

Westgard	تفسیر	Free T3	آزمایش	SIEMENS Centaur CP	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	7.48	SD	0.11	Exp. Date	1402/9/12	Lot. No	36-02-026
Next Month	Mean	7.44	SD	0.1	CV	1.37			

