

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	2.65	
2 - 1	2.72	
4 - 1	2.85	
5 - 1	2.81	
6 - 1	2.89	
7 - 1	2.95	
8 - 1	2.67	
9 - 1	2.69	
11 - 1	2.75	
12 - 1	2.84	
13 - 1	2.85	
14 - 1	2.78	
16 - 1	2.89	
18 - 1	2.91	
19 - 1	2.82	
20 - 1	2.79	
21 - 1	2.84	
22 - 1	2.92	
26 - 1	2.9	
27 - 1	2.86	
28 - 1	2.79	
29 - 1	2.74	
30 - 1	2.81	

Westgard	تفسیر	Free T4	آزمایش	2000 ایمولایت XPI-	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	2.82	SD	0.09	Exp. Date	1403/7/12	Lot. No	36-02-028
Next Month	Mean	2.81	SD	0.08	CV	2.91			

