

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
3 - 1	72.2	1:2S
4 - 1	70.2	
5 - 1	71.4	
6 - 1	70.9	
8 - 1	71.4	
9 - 1	71.2	4:1S
10 - 1	69	
12 - 1	72	
13 - 1	69.6	
15 - 1	69.4	
16 - 1	70.5	
17 - 1	70.8	
20 - 1	68.4	
22 - 1	68.9	
23 - 1	72.2	
24 - 1	68.4	
25 - 1	69.1	
26 - 1	67	1:2S
27 - 1	71	1:2S 2:2S
29 - 1	68.1	
30 - 1	70.3	

Westgard	تفسیر	Gamma GT	آزمایش	اتوآنالایزر BS-800	دستگاه	7	ماه	1402	سال
Control	biorex	Mean	69.5	SD	1.77	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	70.1	SD	1.47	CV	2.09			

