

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 0.99	
Mean	11.99	Mean	17.47
SD	0.13	SD	0.17
CV	1.11	CV	0.99
Day-Run	Value	Day-Run	Value
4 - 1	12	4 - 1	17.2
5 - 1	12.1	5 - 1	17.6
6 - 1	11.9	6 - 1	17.4
8 - 1	12	8 - 1	17.4
9 - 1	11.9	9 - 1	17.7
10 - 1	11.9	10 - 1	17.5
12 - 1	11.7	12 - 1	17.4
13 - 1	11.9	13 - 1	17.8
13 - 2	12	13 - 2	17.5
15 - 1	12	15 - 1	17.3
16 - 1	11.8	16 - 1	17.4
18 - 1	12	18 - 1	17.6
19 - 1	11.9	19 - 1	17.3
20 - 1	12.1	20 - 1	17.6
22 - 1	12.2	22 - 1	17.8
24 - 1	12	24 - 1	17.3
25 - 1	12	25 - 1	17.3
26 - 1	12.2	26 - 1	17.4
27 - 1	12.2	27 - 1	17.5

Westgard	تفسیر	HB	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	7	ماه	1402	سال
Normal Control	Labex	Mean	12	SD	0.15	Exp. Date	1402/6/30	Lot. No	02-2305n

