

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 0.93	
Mean	39.95	Mean	54.92
SD	0.44	SD	0.51
CV	1.11	CV	0.93
Day-Run	Value	Day-Run	Value
4 - 1	40	4 - 1	55.2
5 - 1	40.3	5 - 1	54.6
6 - 1	39.4	6 - 1	55.5
8 - 1	40	8 - 1	55.2
9 - 1	40.6	9 - 1	55.2
10 - 1	40	10 - 1	55.3
12 - 1	40.8	12 - 1	54.5
13 - 1	40	13 - 1	55.6
15 - 1	40	15 - 1	55
16 - 1	39.5	16 - 1	54.2
19 - 1	39.2	19 - 1	54.4
19 - 2	40.1	19 - 2	54.6
20 - 1	40.4	20 - 1	56
22 - 1	40.3	22 - 1	55.1
24 - 1	39.5	24 - 1	54.5
25 - 1	39.4	25 - 1	54.5
26 - 1	39.6	26 - 1	55
27 - 1	40	27 - 1	54.2

Westgard	تفسیر	HCT	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	7	ماه	1402	سال
Normal Control	Labex	Mean	39.8	SD	0.5	Exp. Date	1402/6/30	Lot. No	2305n

