

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal	
Mean	39.82
SD	0.5
CV	1.24
Day-Run	Value
1 - 1	39.2
2 - 1	39
4 - 1	40
5 - 1	40.1
6 - 1	39
7 - 1	39.4
8 - 1	40.3
9 - 1	39.5
11 - 1	40.5
12 - 1	40.5
13 - 1	40
18 - 1	40
19 - 1	40.8
20 - 1	39.5
21 - 1	40
22 - 1	40
28 - 1	39.2
26 - 1	39.5
27 - 1	40
28 - 2	39.5
29 - 1	40
29 - 2	40

High - CV : 1.08	
Mean	54.93
SD	0.59
CV	1.08
Day-Run	Value
1 - 1	55
2 - 1	55
4 - 1	55.5
5 - 1	54.3
6 - 1	55.8
7 - 1	54.4
8 - 1	55.7
9 - 1	55
11 - 1	55
12 - 1	55.5
13 - 1	54
18 - 1	55
19 - 1	55
20 - 1	55
21 - 1	55.5
22 - 1	55.5
28 - 1	54.1
26 - 1	55.6
27 - 1	54.1
28 - 2	54.3
29 - 1	54.2

Westgard	تفسیر	HCT	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	6	ماه	1402	سال
Normal Control	Labex	Mean	39.76	SD	0.5	Exp. Date	1402/6/30	Lot. No	2305n

