

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
3 - 1	65.6	
4 - 1	66.5	
5 - 1	66	
6 - 1	68.1	
8 - 1	65.1	
9 - 1	66	
10 - 1	68.4	
12 - 1	66	
13 - 1	67.2	
15 - 1	64	1:2S
16 - 1	66	
17 - 1	65.5	
20 - 1	65.2	
22 - 1	68	
23 - 1	69	1:2S
24 - 1	65	
25 - 1	66	
26 - 1	67.1	
27 - 1	65.8	
29 - 1	63.8	
30 - 1	68.1	

Westgard	تفسیر	HDL	آزمایش	BS-800	دستگاه	7	ماه	1402	سال
Control	biorex	Mean	66	SD	1.39	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	66.3	SD	1.42	CV	2.14			

