

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	68	
2 - 1	68	
4 - 1	69.7	
5 - 1	67.2	
6 - 1	67	
7 - 1	68.1	
8 - 1	67.6	
11 - 1	67.4	
12 - 1	66	
13 - 1	69	
14 - 1	66	
18 - 1	67.6	
19 - 1	69.1	
20 - 1	64.8	1:2S
21 - 1	68	
22 - 1	67.1	
26 - 1	66.9	
27 - 1	69.7	
28 - 1	65	
29 - 1	66	

Westgard	تفسیر	HDL	آزمایش	اتوآنالایزر BS-800	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	68	SD	1.5	Exp. Date	1402/9/12	Lot. No	035
Next Month	Mean	67.41	SD	1.39	CV	2.06			

