

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	133	
4 - 1	135.6	
4 - 2	140.3	
5 - 1	138.1	
6 - 1	139.2	
8 - 1	138.8	
9 - 1	143	1:2S
10 - 1	139.9	
12 - 1	132.4	
13 - 1	132.7	
15 - 1	129.9	
16 - 1	132.7	
17 - 1	140	
18 - 1	140.6	
19 - 1	140.2	
22 - 1	132.5	
23 - 1	136.8	
24 - 1	138.8	
25 - 1	141.1	
27 - 1	140.6	
28 - 1	140.2	
30 - 1	133.9	

Westgard	تفسیر	Ig A	آزمایش	BS-480 اتوآنالایزر	دستگاه	7	ماه	1402	سال
Control	biorex	Mean	135.79	SD	3.14	Exp. Date	1403/9/14	Lot. No	641-01-030/2
Next Month	Mean	137.29	SD	3.75	CV	2.73			

