

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	46.2	
3 - 1	50.4	
4 - 1	47.2	
5 - 1	47	
6 - 1	50.4	
8 - 1	57.3	
9 - 1	47.2	
10 - 1	49.7	
12 - 1	48.7	
13 - 1	46.3	
15 - 1	51.5	
16 - 1	46.6	
17 - 1	42.4	1:2S
18 - 1	48.9	
19 - 1	46.7	
20 - 1	48.8	
22 - 1	51.5	
23 - 1	48.5	
24 - 1	57	
25 - 1	52.1	
26 - 1	50	
27 - 1	46.5	
29 - 1	44.3	
30 - 1	57.1	

Westgard	تفسیر	Ig M	آزمایش	BS-480 اتوآنالایزر	دستگاه	7	ماه	1402	سال
Control	biorex	Mean	50.05	SD	3.7	Exp. Date	1403/7/7	Lot. No	641-01-030/2
Next Month	Mean	49.26	SD	3.81	CV	7.73			

