

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	51.2	
2 - 1	52.5	
4 - 1	55	Cusum 3S
6 - 1	48	Cusum 3S
7 - 1	53.4	
8 - 1	48	
9 - 1	46.5	
11 - 1	47.1	
12 - 1	54.1	
13 - 1	47.5	
14 - 1	45.2	
16 - 1	53.1	
18 - 1	53	
19 - 1	44.09	
20 - 1	48	
21 - 1	53	
22 - 1	47.8	
26 - 1	55	
27 - 1	44.6	
28 - 1	52.07	
29 - 1	55.3	
30 - 1	46.7	

Westgard	تفسیر	Ig M	آزمایش	BS-480 اتوآنالایزر	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	47.83	SD	3.78	Exp. Date	1403/9/14	Lot. No	641-01-030/2
Next Month	Mean	50.05	SD	3.7	CV	7.4			

