

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	58.6	
3 - 1	62.1	
4 - 1	66.7	
5 - 1	67.03	
6 - 1	61.9	
8 - 1	66.2	
9 - 1	68	
10 - 1	61	
12 - 1	62.7	
13 - 1	59.7	
15 - 1	64.7	
16 - 1	68.2	
17 - 1	65.4	
18 - 1	63	
19 - 1	62.5	
20 - 1	64.7	
22 - 1	71.2	Cusum 3S
23 - 1	66.8	
24 - 1	61.8	
25 - 1	58.6	
26 - 1	58.9	
27 - 1	56.2	
29 - 1	67.5	

Westgard	تفسیر	Insulin (Base)	آزمایش	SIEMENS Centaur CP	دستگاه	7	ماه	1402	سال
Control	biorex	Mean	63.23	SD	4.9	Exp. Date	1403/6/19	Lot. No	36-02-028
Next Month	Mean	63.63	SD	3.8	CV	5.97			

