

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	232.3	
2 - 1	240	
4 - 1	237.8	
5 - 1	239.7	
6 - 1	231	4:1S
6 - 2	230	
7 - 1	241	
8 - 1	240	
11 - 1	225	
12 - 1	239	
13 - 1	235.2	
14 - 1	235	
18 - 1	230	
19 - 1	235.2	
20 - 1	238	
21 - 1	226	1:2S
22 - 1	234.2	
26 - 1	234.5	
27 - 1	227	
28 - 1	235	
29 - 1	231.3	

Westgard	تفسیر	Iron	آزمایش	BS-800 اتوآنالایزر	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	232.95	SD	4.65	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	234.15	SD	4.79	CV	2.05			

