

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	309.54
SD	8.79
CV	2.84
Day-Run	Value
3 - 1	302
4 - 1	313.4
5 - 1	300
6 - 1	318
8 - 1	319
9 - 1	316
10 - 1	301
12 - 1	320
13 - 1	307
15 - 1	298
16 - 1	312
17 - 1	301
20 - 1	313
22 - 1	300
23 - 1	320
24 - 1	312
25 - 1	294
26 - 1	319
27 - 1	320
29 - 1	300
30 - 1	315

High - CV : 1.57	
Mean	433.5
SD	6.79
CV	1.57
Day-Run	Value
3 - 1	433
4 - 1	437.6
5 - 1	420
6 - 1	441
8 - 1	450
9 - 1	428
10 - 1	430
12 - 1	435
13 - 1	428
15 - 1	431
16 - 1	439
17 - 1	428
20 - 1	434
22 - 1	438
23 - 1	436
24 - 1	438
25 - 1	425
26 - 1	442
27 - 1	428
29 - 1	428
30 - 1	434

Westgard	تفسیر	LDH	آزمایش	اتوآنالایزر BS-800	دستگاه	7	ماه	1402	سال
Normal Control	biorex	Mean	310.28	SD	9.47	Exp. Date	1403/1/17	Lot. No	035

