

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 2.06	
Mean	310.28	Mean	436.43
SD	9.47	SD	8.99
CV	3.05	CV	2.06
Day-Run	Value	Day-Run	Value
1 - 1	310	1 - 1	424
2 - 1	311	2 - 1	427
4 - 1	313	4 - 1	448
5 - 1	301	5 - 1	439
6 - 1	330	6 - 1	423
8 - 1	299	8 - 1	449
11 - 1	308	11 - 1	426
12 - 1	324.4	12 - 1	445.2
13 - 1	308	13 - 1	433
14 - 1	316	14 - 1	425
18 - 1	307	18 - 1	436
19 - 1	295	19 - 1	430
20 - 1	315	20 - 1	451
21 - 1	298	21 - 1	436
22 - 1	314	22 - 1	439
26 - 1	300	26 - 1	440
27 - 1	323	27 - 1	434
28 - 1	306	28 - 1	447
29 - 1	317	29 - 1	440

Westgard	تفسیر	LDH	آزمایش	اتوآنالایزر BS-800	دستگاه	6	ماه	1402	سال
Normal Control	biorex	Mean	313.05	SD	9.45	Exp. Date	1403/1/17	Lot. No	035

