

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	135.8	
2 - 1	130	1:2S
4 - 1	140	
5 - 1	134.9	
6 - 1	138	
7 - 1	137.8	
8 - 1	138	
11 - 1	133	
12 - 1	139	
13 - 1	138.4	
14 - 1	140	
18 - 1	141	
19 - 1	136.1	
20 - 1	136.6	
21 - 1	139	
22 - 1	138.8	
26 - 1	140	
27 - 1	135	
28 - 1	141	
29 - 1	138.4	

Westgard	تفسیر	LDL	آزمایش	اتوآنالایزر BS-800	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	137.19	SD	2.53	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	137.54	SD	2.78	CV	2.02			

