

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	32.4	
3 - 1	33.46	
4 - 1	33.52	
5 - 1	33.87	
6 - 1	32.46	
7 - 1	32.5	
8 - 1	32.8	
9 - 1	33.5	
10 - 1	31.4	
12 - 1	34.1	
13 - 1	33.9	
15 - 1	32.9	
16 - 1	34.4	
17 - 1	32.8	
18 - 1	31.4	
19 - 1	34.9	
20 - 1	32.6	
22 - 1	32.8	
23 - 1	35.2	1:2S
24 - 1	33.2	1:3S
25 - 1	33.4	
26 - 1	32.8	
27 - 1	31.2	1:3S
29 - 1	32.6	
30 - 1	31.2	1:2S

Westgard	تفسیر	LH	آزمایش	2000 ایمولایت XPI-	دستگاه	7	ماه	1402	سال
Control	Seronorm	Mean	32.64	SD	1.27	Exp. Date	1403/5/31	Lot. No	2011909
Next Month	Mean	33	SD	1.23	CV	3.74			

