

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	43.2	
3 - 1	44.1	
4 - 1	43.5	
5 - 1	43.2	
6 - 1	43.5	
8 - 1	42.73	
9 - 1	43.5	
9 - 2	43	
10 - 1	43.5	
12 - 1	44	
13 - 1	42.8	
15 - 1	43.4	
16 - 1	43.9	
17 - 1	43.8	
18 - 2	43.2	
19 - 1	43.3	
20 - 1	44.2	1:2S
22 - 1	43.3	
23 - 1	43.8	
24 - 1	42.9	
25 - 1	44.1	
26 - 1	43.5	
27 - 1	42.4	1:2S
29 - 1	43.3	
30 - 1	43.5	

Westgard	تفسیر	Lactate(Lactic Acid)	آزمایش	BS-480 اتوآنالایزر	دستگاه	7	ماه	1402	سال
Control	biorex	Mean	43.33	SD	0.42	Exp. Date	1402/12/14	Lot. No	
Next Month	Mean	43.43	SD	0.45	CV	1.05			

