

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	44.2	
2 - 1	43.5	
4 - 1	43.6	
5 - 1	43.6	
6 - 1	43.4	Cusum 3S
8 - 1	43.3	
9 - 1	43.9	
11 - 1	42.8	
12 - 1	43.2	
13 - 1	43.5	
14 - 1	42.6	1:3S
16 - 1	43.3	
18 - 1	42.7	
19 - 1	43.1	
20 - 1	43.3	
21 - 1	43.9	
22 - 1	42.9	
26 - 1	43.1	
30 - 1	43.4	

Westgard	تفسیر	Lactate(Lactic Acid)	آزمایش	BS-480 اتوآنالایزر	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	43.22	SD	0.69	Exp. Date	1402/12/14	Lot. No	
Next Month	Mean	43.33	SD	0.42	CV	0.96			

