

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
3 - 1	67.4	
4 - 1	63.4	1:2S
5 - 1	68.2	
6 - 1	64.2	
8 - 1	68.2	
9 - 1	69.1	
10 - 1	64.2	
12 - 1	64.7	
13 - 1	64.7	
15 - 1	67.9	
16 - 1	69.7	
17 - 1	65.1	
20 - 1	69.7	
22 - 1	64.2	
23 - 1	66.4	
24 - 1	69.9	
25 - 1	66.5	
26 - 1	68	
27 - 1	69.7	
29 - 1	69.1	
30 - 1	65	

Westgard	تفسیر	Lipase	آزمایش	اتوآنالایزر BS-800	دستگاه	7	ماه	1402	سال
Control	biorex	Mean	66.2	SD	1.85	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	66.92	SD	2.23	CV	3.33			

