

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 0.78	
Mean	30.18	Mean	31.83
SD	0.87	SD	0.25
CV	2.89	CV	0.78
Day-Run	Value	Day-Run	Value
4 - 1	31.5	4 - 1	31.5
5 - 1	29.9	5 - 1	31.9
6 - 1	30.5	6 - 1	32
8 - 1	29.3	8 - 1	31.5
9 - 1	29.5	9 - 1	32
10 - 1	30.6	10 - 1	32.2
12 - 1	31.5	12 - 1	32
13 - 1	29	13 - 1	32
15 - 1	29.8	15 - 1	31.7
16 - 1	30	16 - 1	31.5
18 - 1	31	18 - 1	32
19 - 1	29.3	19 - 1	32.1
20 - 1	29.1	20 - 1	31.4
22 - 1	29.1	22 - 1	31.9
24 - 1	29.9	24 - 1	32
25 - 1	31	25 - 1	31.6
26 - 1	31.2	26 - 1	32
27 - 1	31	27 - 1	31.6

Westgard	تفسیر	MCHC	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	7	ماه	1402	سال
Normal Control	Labex	Mean	30.32	SD	0.65	Exp. Date	1402/6/30	Lot. No	02-2305n

