

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal	
Mean	92.03
SD	0.7
CV	0.76
Day-Run	Value
4 - 1	92
5 - 1	92.5
6 - 1	91.3
8 - 1	92.8
9 - 1	92.3
10 - 1	91
12 - 1	92.5
13 - 1	93
15 - 1	92.3
16 - 1	91.7
18 - 1	91.3
19 - 1	93
20 - 1	92.5
22 - 1	90.7
24 - 1	92
25 - 1	91.1
26 - 1	92.2
27 - 1	92.4

High - CV : 0.61	
Mean	99.41
SD	0.6
CV	0.61
Day-Run	Value
4 - 1	98
5 - 1	99.5
6 - 1	99.2
8 - 1	100
9 - 1	99
10 - 1	100.1
12 - 1	99
13 - 1	99.5
15 - 1	98.3
16 - 1	99.1
18 - 1	100.1
19 - 1	99.2
20 - 1	99.3
22 - 1	100
24 - 1	99.4
25 - 1	100
26 - 1	99.6
27 - 1	100

Westgard	تفسیر	MCV	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	7	ماه	1402	سال
Normal Control	Labex	Mean	91.85	SD	0.6	Exp. Date	1402/6/30	Lot. No	02-2305n

