

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

| Normal | |
|---------|-------|
| Mean | 91.85 |
| SD | 0.6 |
| CV | 0.65 |
| Day-Run | Value |
| 1 - 1 | 91 |
| 2 - 1 | 92.2 |
| 4 - 1 | 92 |
| 5 - 1 | 91.7 |
| 6 - 1 | 91.2 |
| 7 - 1 | 91.2 |
| 8 - 1 | 92.3 |
| 9 - 1 | 91 |
| 11 - 1 | 92 |
| 12 - 1 | 92 |
| 13 - 1 | 92.5 |
| 14 - 1 | 92.3 |
| 18 - 1 | 90.9 |
| 20 - 1 | 92.5 |
| 20 - 2 | 93 |
| 21 - 1 | 92.3 |
| 22 - 1 | 92 |
| 28 - 1 | 92.3 |
| 26 - 1 | 91 |
| 27 - 1 | 91.3 |
| 28 - 2 | 92 |
| 29 - 1 | 92 |

| High - CV : 0.61 | |
|------------------|-------|
| Mean | 99.15 |
| SD | 0.6 |
| CV | 0.61 |
| Day-Run | Value |
| 1 - 1 | 98 |
| 2 - 1 | 99 |
| 4 - 1 | 99.5 |
| 5 - 1 | 99 |
| 6 - 1 | 99.5 |
| 7 - 1 | 99 |
| 8 - 1 | 100 |
| 9 - 1 | 99 |
| 11 - 1 | 99.5 |
| 12 - 1 | 99 |
| 13 - 1 | 99.5 |
| 14 - 1 | 98.5 |
| 18 - 1 | 99 |
| 20 - 1 | 100 |
| 20 - 2 | 99.2 |
| 21 - 1 | 99.5 |
| 22 - 1 | 99 |
| 28 - 1 | 100 |
| 26 - 1 | 98.5 |
| 27 - 1 | 98 |
| 28 - 2 | 100 |
| 29 - 1 | 98.5 |

| | | | | | | | | | |
|----------------|-------|------|--------|-----------------------|--------|-----------|-----------|---------|----------|
| Westgard | تفسیر | MCV | آزمایش | NIHON KOHDEN سل کانتر | دستگاه | 6 | ماه | 1402 | سال |
| Normal Control | Labex | Mean | 91.65 | SD | 0.64 | Exp. Date | 1402/6/30 | Lot. No | 02-2305n |

