

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	2.29
SD	0.16
CV	7.09
Day-Run	Value
3 - 1	2
4 - 1	2
5 - 1	2.3
6 - 1	2.3
8 - 1	2.4
9 - 1	2.3
10 - 1	2.3
12 - 1	2.5
13 - 1	2.1
15 - 1	2.5
16 - 1	2.4
17 - 1	2.2
20 - 1	2.4
22 - 1	2
23 - 1	2.5
24 - 1	2.4
25 - 1	2.4
26 - 1	2.1
27 - 1	2.3
29 - 1	2.3
30 - 1	2.3

High - CV : 5.09	
Mean	3.76
SD	0.19
CV	5.09
Day-Run	Value
3 - 1	3.8
4 - 1	3.8
5 - 1	3.9
6 - 1	3.5
8 - 1	4
9 - 1	3.7
10 - 1	3.8
12 - 1	3.9
13 - 1	3.7
15 - 1	4
16 - 1	3.4
17 - 1	4
20 - 1	3.6
22 - 1	4
23 - 1	3.7
24 - 1	3.5
25 - 1	3.8
26 - 1	3.9
27 - 1	3.7
29 - 1	3.4
30 - 1	3.8

Westgard	تفسیر	Magnesium(Mg)	آزمایش	اتوآنالایزر BS-800	دستگاه	7	ماه	1402	سال
Normal Control	biorex	Mean	2.22	SD	0.14	Exp. Date	1403/1/17	Lot. No	035

