

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 2.31	
Mean	226.76	Mean	549
SD	3.49	SD	12.7
CV	1.54	CV	2.31
Day-Run	Value	Day-Run	Value
4 - 1	225	4 - 1	537
5 - 1	220	5 - 1	544
6 - 1	229	6 - 1	560
8 - 1	221	8 - 1	533
9 - 1	225	9 - 1	545
10 - 1	227.2	10 - 1	573
13 - 1	230	13 - 1	535
15 - 1	232	15 - 1	550
16 - 1	227.8	16 - 1	561
18 - 1	225	18 - 1	554
19 - 1	230	19 - 1	554
20 - 1	224	20 - 1	533
22 - 1	229	22 - 1	534
24 - 1	231	24 - 1	558
25 - 1	224	25 - 1	542
26 - 1	230	26 - 1	570
27 - 1	225	27 - 1	550

Westgard	تفسیر	PLT	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	7	ماه	1402	سال
Normal Control	Labex	Mean	224.85	SD	3.84	Exp. Date	1402/6/30	Lot. No	02-2305n

