

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 1.76	
Mean	224.85	Mean	551.55
SD	3.84	SD	9.71
CV	1.71	CV	1.76
Day-Run	Value	Day-Run	Value
1 - 1	230	1 - 1	540
2 - 1	225	2 - 1	550
4 - 1	221	4 - 1	560
5 - 1	225	5 - 1	550
6 - 1	222	6 - 1	545
7 - 1	225	7 - 1	560
8 - 1	228	8 - 1	543
9 - 1	220	9 - 1	556
11 - 1	221	11 - 1	540
12 - 1	225	12 - 1	560
13 - 1	225	13 - 1	535
14 - 1	230	14 - 1	540
18 - 1	220	18 - 1	570
19 - 1	229	19 - 1	560
20 - 1	225	20 - 1	550
22 - 1	217	22 - 1	559
28 - 1	230	28 - 1	560
27 - 1	230	27 - 1	540
28 - 2	224	28 - 2	552
29 - 1	225	29 - 1	561

Westgard	تفسیر	PLT	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	6	ماه	1402	سال
Normal Control	Labex	Mean	223.39	SD	3.76	Exp. Date	1402/6/30	Lot. No	02-2305n

