

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	3.24	
3 - 1	3.52	
4 - 1	3.48	
5 - 1	3.62	
6 - 1	3.21	
7 - 1	3.4	
8 - 1	3.7	
9 - 1	3.4	
10 - 1	3.3	
12 - 1	3.2	
13 - 1	3.7	
15 - 1	3.4	
16 - 1	3.4	
17 - 1	2.99	
18 - 1	3.1	
19 - 1	3.6	
20 - 1	3.43	
22 - 1	3.12	
23 - 1	3.16	
24 - 1	3.28	
25 - 1	3.45	
26 - 1	3.56	
27 - 1	3.45	
29 - 1	2.98	
30 - 1	3.08	

Westgard	تفسیر	PSA	آزمایش	2000 ایمولایت XPI-	دستگاه	7	ماه	1402	سال
Control	biorex	Mean	3.35	SD	0.19	Exp. Date	1402/7/1	Lot. No	365-02-028
Next Month	Mean	3.35	SD	0.21	CV	6.27			

