

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 3.47	
Mean	3.51	Mean	8.39
SD	0.09	SD	0.29
CV	2.69	CV	3.47
Day-Run	Value	Day-Run	Value
1 - 1	3.5	1 - 1	8.6
2 - 1	3.4	2 - 1	8.4
4 - 1	3.6	4 - 1	8.9
5 - 1	3.5	5 - 1	8.6
6 - 1	3.6	6 - 1	8
7 - 1	3.6	7 - 1	8.4
8 - 1	3.6	8 - 1	8
11 - 1	3.4	11 - 1	8.5
12 - 1	3.4	12 - 1	8.1
13 - 1	3.4	13 - 1	8.3
14 - 1	3.4	14 - 1	8.6
18 - 1	3.5	18 - 1	8
19 - 1	3.5	19 - 1	8.4
20 - 1	3.4	20 - 1	8.1
21 - 1	3.4	21 - 1	8.4
22 - 1	3.5	22 - 1	8.6
26 - 1	3.6	26 - 1	9
27 - 1	3.5	27 - 1	8
28 - 1	3.7	28 - 1	8.4
29 - 1	3.6	29 - 1	8.4

Westgard	تفسیر	Phosphorus (P)	آزمایش	اتوآنالایزر BS-800	دستگاه	6	ماه	1402	سال
Normal Control	biorex	Mean	3.45	SD	0.11	Exp. Date	1403/1/17	Lot. No	035

