

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	18.1	
3 - 1	18.72	
4 - 1	18.96	
5 - 1	18.42	
6 - 1	17.94	
7 - 1	17.7	
8 - 1	18.4	
9 - 1	18.2	
10 - 1	19.2	
12 - 1	19.1	
13 - 1	18.7	
15 - 1	18.5	
16 - 1	18.5	
17 - 1	19.6	
18 - 1	18.7	
19 - 1	18.9	
20 - 1	19.2	
22 - 1	17.5	
23 - 1	17.9	1:3S
24 - 1	18.2	1:2S
25 - 1	18.2	
26 - 1	19.3	
27 - 1	18.5	
29 - 1	16.8	1:2S
30 - 1	17.3	1:2S 2:2S

Westgard	تفسیر	Prolactin	آزمایش	2000 ایمولایت XPI-	دستگاه	7	ماه	1402	سال
Control	biorex	Mean	18.64	SD	0.62	Exp. Date	1403/2/10	Lot. No	365-02-028
Next Month	Mean	18.42	SD	0.67	CV	3.61			

