

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 1.42	
Mean	4.35	Mean	5.51
SD	0.04	SD	0.08
CV	1	CV	1.42
Day-Run	Value	Day-Run	Value
4 - 1	4.31	4 - 1	5.6
5 - 1	4.32	5 - 1	5.52
6 - 1	4.35	6 - 1	5.56
8 - 1	4.3	8 - 1	5.48
9 - 1	4.43	9 - 1	5.62
10 - 1	4.38	10 - 1	5.5
12 - 1	4.29	12 - 1	5.49
13 - 1	4.35	13 - 1	5.6
13 - 2	4.38	13 - 2	5.4
15 - 1	4.35	15 - 1	5.5
16 - 1	4.29	16 - 1	5.4
18 - 1	4.42	18 - 1	5.58
19 - 1	4.32	19 - 1	5.51
20 - 1	4.39	20 - 1	5.62
22 - 1	4.33	22 - 1	5.59
24 - 1	4.4	24 - 1	5.41
25 - 1	4.3	25 - 1	5.38
26 - 1	4.38	26 - 1	5.5
27 - 1	4.35	27 - 1	5.45

Westgard	تفسیر	RBC	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	7	ماه	1402	سال
Normal Control	Labex	Mean	4.36	SD	0.05	Exp. Date	1402/6/30	Lot. No	02-2305n

