

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	35	
2 - 1	38.2	
5 - 1	34.4	
6 - 1	36.7	
7 - 1	35.4	
8 - 1	37.3	
9 - 1	36.5	1:3S
11 - 2	35.8	
12 - 1	37.2	
13 - 1	38.5	
14 - 1	33.5	
16 - 1	34.5	
18 - 1	38.2	1:3S
19 - 1	36.9	
20 - 1	37	
21 - 1	36.5	
22 - 1	38.1	
27 - 1	37.4	
28 - 1	36.1	
29 - 1	36.6	
30 - 1	37.4	

Westgard	تفسیر	RF	آزمایش	BS-480 اتوآنالایزر	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	35.92	SD	1.46	Exp. Date	1405/9/14	Lot. No	645-03-018/2
Next Month	Mean	36.53	SD	1.36	CV	3.73			

