

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	145	
2 - 1	146.5	
4 - 1	148	
5 - 1	146.5	
6 - 1	147	
7 - 1	146.3	
8 - 1	147.2	
9 - 1	145	
11 - 1	146.2	
12 - 1	146.8	
12 - 2	146.1	
14 - 1	146.3	
16 - 1	146.9	
18 - 1	147	
19 - 1	147.5	
20 - 1	146.3	
21 - 1	146.8	
22 - 1	146.1	
26 - 1	147.5	
27 - 1	146.1	

Westgard	تفسیر	Sodium ( Na )	آزمایش	AUDICOM	دستگاه	6	ماه	1402	سال
Control	زیست شیمی	Mean	146.8	SD	1.72	Exp. Date	1404/1/14	Lot. No	761LOT17
Next Month	Mean	146.56	SD	0.75	CV	0.51			

