

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	10.2	
2 - 1	10.9	
4 - 1	12.05	
5 - 1	11.1	
7 - 1	11	
8 - 1	11.36	
9 - 1	11.5	
11 - 1	12	
12 - 1	10.5	
13 - 1	11.62	
14 - 1	11.1	
16 - 1	10.58	
18 - 1	11	
19 - 1	11.52	
20 - 1	11.7	1:3S
21 - 1	10.96	
22 - 1	10.63	
26 - 1	11.42	
27 - 1	12	
28 - 1	11.35	
29 - 1	11.2	
30 - 1	10.65	

Westgard	تفسیر	T4	آزمایش	SIEMENS Centaur CP	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	11.07	SD	0.53	Exp. Date	1403/6/19	Lot. No	36-02-028
Next Month	Mean	11.2	SD	0.51	CV	4.56			

