

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	315	
2 - 1	315.3	
4 - 1	324	
5 - 1	316	
6 - 1	328	
7 - 1	316.2	
8 - 1	322.9	
11 - 1	324	
12 - 1	307	
13 - 1	315.4	
14 - 1	323.3	
18 - 1	313.1	
19 - 1	316	
20 - 1	314	
21 - 1	285	1:2S
22 - 1	317.5	
26 - 1	314.4	
27 - 1	320	
28 - 1	310	
29 - 1	308.5	

Westgard	تفسیر	TIBC	آزمایش	اتوآنالایزر BS-800	دستگاه	6	ماه	1402	سال
Control	biorex	Mean	315.84	SD	10.39	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	315.28	SD	8.98	CV	2.85			

