

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	3.35	
3 - 1	3.64	
4 - 1	3.8	
5 - 1	3.97	
6 - 1	3.73	
8 - 1	3.84	
9 - 1	3.95	
10 - 1	4	
12 - 1	3.37	
13 - 1	3.45	
15 - 1	3.89	
16 - 1	3.63	
17 - 1	3.59	
18 - 1	3.9	
19 - 1	3.65	
20 - 1	3.59	
22 - 1	3.84	
23 - 1	3.89	
24 - 1	3.53	1:3S
26 - 1	3.45	
27 - 1	3.6	
29 - 1	3.93	

Westgard	تفسیر	TSH	آزمایش	SIEMENS Centaur CP	دستگاه	7	ماه	1402	سال
Control	biorex	Mean	3.67	SD	0.2	Exp. Date	1403/6/19	Lot. No	36-02-028
Next Month	Mean	3.71	SD	0.2	CV	5.49			

