

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 1.54	
Mean	171.19	Mean	233.05
SD	1.72	SD	3.6
CV	1.01	CV	1.54
Day-Run	Value	Day-Run	Value
3 - 1	173	3 - 1	226
4 - 1	167	4 - 1	233
5 - 1	170	5 - 1	235
6 - 1	170	6 - 1	230
8 - 1	170	8 - 1	230
9 - 1	174	9 - 1	240
10 - 1	172	10 - 1	236
12 - 1	174	12 - 1	230
13 - 1	170	13 - 1	232.1
15 - 1	170	15 - 1	240
16 - 1	170	16 - 1	232
17 - 1	171	17 - 1	235
20 - 1	173	20 - 1	230
22 - 1	170	22 - 1	233
23 - 1	171	23 - 1	235
24 - 1	173	24 - 1	235
25 - 1	171	25 - 1	231
26 - 1	170	26 - 1	237
27 - 1	171	27 - 1	230
29 - 1	173	29 - 1	235
30 - 1	172	30 - 1	229

Westgard	تفسیر	Triglycerides	آزمایش	اتوآنالایزر BS-800	دستگاه	7	ماه	1402	سال
Normal Control	biorex	Mean	170.53	SD	2.34	Exp. Date	1403/1/17	Lot. No	035

