

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.47	
Mean	34.36	Mean	132.91
SD	1.05	SD	1.95
CV	3.05	CV	1.47
Day-Run	Value	Day-Run	Value
1 - 1	34.7	1 - 1	132.7
2 - 1	35.2	2 - 1	132.5
4 - 1	34.2	4 - 1	133.3
5 - 1	34.3	5 - 1	132.8
6 - 1	35.9	6 - 1	135.2
7 - 1	33.6	7 - 1	130.6
8 - 1	35.8	8 - 1	133
11 - 1	34.5	11 - 1	130
12 - 1	33	12 - 1	132.6
13 - 1	33.6	13 - 1	135.2
14 - 1	36	14 - 1	130
18 - 1	34.4	18 - 1	135.6
19 - 1	35.3	19 - 1	133.6
20 - 1	35.2	20 - 1	130.3
21 - 1	33.2	21 - 1	134.9
22 - 1	33.5	22 - 1	134.3
26 - 1	33.8	26 - 1	133
27 - 1	35	27 - 1	132.4
28 - 1	32.1	28 - 1	130
29 - 1	33.8	29 - 1	136.1

Westgard	تفسیر	Urea	آزمایش	BS-800 اتوآنالایزر	دستگاه	6	ماه	1402	سال
Normal Control	biorex	Mean	34.01	SD	1.02	Exp. Date	1403/1/17	Lot. No	035

