

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 2.52	
Mean	4.16	Mean	8.63
SD	0.1	SD	0.22
CV	2.34	CV	2.52
Day-Run	Value	Day-Run	Value
3 - 1	4.2	3 - 1	8.6
4 - 1	4.2	4 - 1	8.4
5 - 1	4.1	5 - 1	8.6
6 - 1	4.1	6 - 1	8.5
8 - 1	4.1	8 - 1	8.5
9 - 1	4.2	9 - 1	8.5
10 - 1	4.2	10 - 1	8.7
12 - 1	4.3	12 - 1	8.6
13 - 1	4.2	13 - 1	8.2
15 - 1	4.1	15 - 1	8.7
16 - 1	4.2	16 - 1	8.7
17 - 1	4.1	17 - 1	8.5
20 - 1	4	20 - 1	8.5
22 - 1	4	22 - 1	8.8
23 - 1	4.2	23 - 1	8.5
24 - 1	4.1	24 - 1	8.7
25 - 1	4	25 - 1	8.9
26 - 1	4.2	26 - 1	8.7
27 - 1	4.3	27 - 1	9
29 - 1	4.3	29 - 1	8.5
30 - 1	4.3	30 - 1	9.2

Westgard	تفسیر	Uric Acid	آزمایش	اتوانالایزر BS-800	دستگاه	7	ماه	1402	سال
Normal Control	biorex	Mean	4.2	SD	0.1	Exp. Date	1403/1/17	Lot. No	035

