

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 2.86	
Mean	4.2	Mean	8.82
SD	0.1	SD	0.25
CV	2.32	CV	2.86
Day-Run	Value	Day-Run	Value
1 - 1	4.3	1 - 1	9
2 - 1	4.2	2 - 1	8.7
4 - 1	4.2	4 - 1	8.8
5 - 1	4.2	5 - 1	9.2
6 - 1	4.2	6 - 1	8.6
7 - 1	4.4	7 - 1	8.7
8 - 1	4.2	8 - 1	8.5
11 - 1	4	11 - 1	8.8
12 - 1	4.2	12 - 1	8.6
13 - 1	4.3	13 - 1	9
14 - 1	4.1	14 - 1	8.5
18 - 1	4.2	18 - 1	8.8
19 - 1	4.2	19 - 1	8.8
20 - 1	4.2	20 - 1	8.7
21 - 1	4.1	21 - 1	9
22 - 1	4.3	22 - 1	9
26 - 1	4.2	26 - 1	9
27 - 1	4.2	27 - 1	8.9
28 - 1	4.3	28 - 1	8.6
29 - 1	4	29 - 1	9.1
29 - 2	4.3	29 - 2	8.5

Westgard	تفسیر	Uric Acid	آزمایش	BS-800 اتوانالایزر	دستگاه	6	ماه	1402	سال
Normal Control	biorex	Mean	4.2	SD	0.15	Exp. Date	1403/1/17	Lot. No	035

