

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
3 - 1	28.5	
4 - 1	29.5	
5 - 1	29.7	
6 - 1	28.7	
8 - 1	28.9	
9 - 1	28.9	
10 - 1	29.5	
12 - 1	29.4	
13 - 1	28.8	
15 - 1	30	
16 - 1	29.7	
17 - 1	29.5	
18 - 1	28.7	
19 - 1	29.7	
20 - 1	29.5	
22 - 1	29	
23 - 1	28.8	
24 - 1	29.1	
25 - 1	29.5	
26 - 1	29	
27 - 1	28.8	
29 - 1	29.3	
30 - 1	29.6	

Westgard	تفسیر	V D3	آزمایش	HPLC Agilent 1100 -1	دستگاه	7	ماه	1402	سال
Control	آزمایشگاه	Mean	29.32	SD	0.68	Exp. Date	1402/7/30	Lot. No	
Next Month	Mean	29.22	SD	0.41	CV	1.42			

